



Peel Down to Open



DO NOT USE IF PACKAGE IS OPENED OR DAMAGED



perfecto™

STERILE SURGICAL GLOVES
POWDERED

SIZE

5.5

Quality as per ASTM D 3577/ EN 455-1,2 & 3
Powdered using absorbable corn starch U.S.P.
Anatomically shaped-Curved fingers-for user comfort
Micro textured surface for better grip
Do not re sterilize



ONE PAIR

Batch No. :

Mfg. Date :

Exp. Date :

Sterile. Date :



8 9 0 8 0 0 0 8 4 3 0 8 0
ISO 9001:2015 ISO 13485:2016



EC Representative :
AMSTERMED B.V.
Saturnusstraat 46-62,
Unit 022, 2132 HB Hoofddorp,
The Netherlands.
Phone : +31 23 5656337.
Web: www.amstermed.nl

This product contains natural rubber latex which may cause allergic reaction including anaphylactic responses in some individuals. In case of allergy discontinue use and consult a physician. Sterile unless package is damaged or opened. This product is for single use only. After donning, remove powder by wiping gloves thoroughly with a sterile wet sponge or other effective methods. Store in cool dry place away from direct sunlight and ozone. Risk of reuse : May cause infection, allergic reaction and poor barrier protection.

Manufactured By :
Beta Healthcare Products (P) Ltd
Plot No - 21B
Cochin Special Economic Zone
Kakkanad, Cochin - 682 037
+91-484-2413340
boney@betahealthcare.com
www.betahealthcare.com

MADE IN INDIA



Recyclable Material

BL/S/25

Rev. No :01