



TECHNICAL FEATURES

- ▶ Anatomically shaped curved fingers
- ▶ Beaded cuff end
- ▶ Micro textured for better grip
- ▶ Powdered using absorbable Corn Starch USP
- ▶ Sterilized with Ethylene Oxide

WARNINGS

- ▶ This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses in some individuals
- ▶ Sterility guaranteed if package is undamaged or unopened
- ▶ Do not use if the package is damaged
- ▶ After donning, remove the powder by wiping gloves thoroughly with a sterile wet sponge or other effective methods
- ▶ Single use
- ▶ Store in a cool, dry place, away from direct sunlight & ozone.
- ▶ Risk of reuse : May cause infection, allergic reaction and poor barrier protection.

POWDER FREE

pristeen[®]

STERILE SURGICAL GLOVES

6

ISO 9001:2015

ISO 13485:2016

M.R.P. : Rs.
(Inclusive of all taxes)

Manufactured By :

Beta Healthcare Products (P) Ltd
Plot No - 21B, Cochin Special Economic Zone
Kakkanad, Cochin - 682 037, Kerala, INDIA

+91 484-2413340

boney@betahealthcare.com

www.betahealthcare.com

EC REP

Zuverlässige GrößtH
Fürther Straße 27
90429 Nürnberg
Germany

CE 2460 SIZE 6-1^{PAIR}



LOT

LATEX



STERILE EO

+91-484 2413340

boney@betahealthcare.com



www.betahealthcare.com